

Free text editor

The Free Text Editor is used in various parts of the MediSign application and offers the ability to insert free text in **basic formatting**.

Prioritizing **secure operation** in an internet environment, it allows input and formatting for the following:

- Paragraph
- Heading
- List
- Tables

Text with basic formatting can also be imported directly from Word (or another Word Processor) with Copy and Paste.

Free text

Patient overview - Important notes

Joe Bailey

Father name: Ralph

Age: 65 years SSN: 818-60-9519 Code: 32

Frequent episodes of chest pain and discomfort. Coronary angiography with normal findings.

T ▼

B

i

☞

↵

U

X'

X

A

Save





Clear





Exit


[Need help?](#)

Tools

By selecting a piece of text, you have the following tools available:








Moving elements


You can easily move up or down a paragraph (or any other element) or delete it.



Tables

Please read [Free Text Editor - Tables](#)

Example:

 **Important notes**

The patient comes in for a routine checkup for coronary artery disease.

Heading1	Heading2	Heading3
text A1 with formatting	text B1	text C1
text A2	text B2	text C2
text A3	text B3	text C3

Related topics

- [Free Text Editor - Tables](#)
- [History notes templates](#)

From:
<https://docs.medesign.com/> - **docs.medisign.com**

Permanent link:
<https://docs.medisign.com/manual/free-text-editor>

Last update: **2022/09/19 11:21**

