

Uploads

Files uploading is a common and very useful process.

Simple steps:

- Click the area in the dotted box (**point 1**) and then select files (or drag the files you want to upload into it)
- In each case, there are special restrictions (file type, size, etc.) for which you are informed in detail (**point 2**)
- When you are ready, press the "Send" button (**point 3**)

Upload files

Patient overview

Joe Bailey

Father name: Ralph

Age: 65 years SSN: 818-60-9519 Code: 32

1

Drop files here to upload

2

Up to 5 files per upload. Allowed file types: jpeg,jpg,png,pdf. Maximum file size: 10MB.
Remember to check the files for viruses before uploading them.

3

Sent

Cancel

[Add file](#) [Remove all](#)

[Need help?](#)

From:
<https://docs.medesign.com/> - docs.medesign.com

Permanent link:
<https://docs.medesign.com/manual/uploads>

Last update: 2021/08/14 08:52

